# TECHNICAL NOTE ABSOLUTE QUANTUM YIELD OF UV-TO-NIR EMITTING SAMPLES

TN\_P19 v.1; 15 Sept. 15, Georgios Arnaoutakis

### INTRODUCTION

The photoluminescence quantum yield (QY) or quantum efficiency  $\eta$  is one of the core parameters in the characterisation of luminescent materials<sup>1,2</sup>. The absolute QY, defined as the ratio of the photons emitted by the material over the absorbed photons, has widely replaced the relative QY method, requiring a fluorescent standard of known QY. It also constitutes a simpler method relying on minimum alignment and calibration compared to thermal



FS5 Spectroflurometer

lens/beam deflection and photoacoustic methods<sup>3</sup>. Instead, a fluorescence spectrometer equipped with an integrating sphere is employed in the absolute method.

## METHODS AND MATERIALS

Liquid and solid materials were measured spanning excitation and emission wavelengths from the UV to the NIR. LED phosphor powders and thin films doped with lanthanides such as Cerium (Ce<sup>3+</sup>), Terbium (Tb<sup>3+</sup>) and Europium (Eu<sup>2+</sup>) were measured in standard PTFE trays. In the NIR range, semiconductor nanocrystal lead sulfide (PbS) quantum dots (QD-NIR-1V, Ocean Optics) in toluene were dispensed into 10 mm path-length quartz cuvettes.

Excitation and emission spectra were measured in a fluorescence spectrometer FS5-NIR equipped with an integrating sphere module (SC-30) and a single photon PMT detector (Hamamatsu, R2658P) extending the emission detection range to 1010 nm. For liquid samples a measurement over the scattering  $L_{sam}$  and emission  $E_{sam}$  of the sample, followed by a measurement of the solvent, also called reference or blank,  $L_{ref}$  and  $E_{ref'}$  as:

$$QY = \frac{E_{sam} - E_{ref}}{L_{ref} - L_{sam}}$$

For solid powder and thin film samples three measurements were performed for the calculation of the quantum yield: i) A measurement of the scattering without the sample,  $L_{ref'}$  ii) a measurement of the emission of the sample under direct excitation,  $E_{samdir'}$  and iii) a measurement of the scattering with the sample under indirect excitation,  $E_{samid'}$  as<sup>4</sup>:



**EDINBURGH** 

INSTRUMENTS

The calculation of the QY was performed in the wizard of the instrument's operating software, Fluoracle.

## RESULTS - DISCUSSION

Figures 1 and 2 display the scattering and emission spectra of Ce<sup>3+</sup> and Ce<sup>3+</sup>/Tb<sup>3+</sup>-doped phosphor, respectively.

In addition to the intense green emission, the QY increases from 85.1% for the Ce<sup>3+</sup>-doped phosphor to 92% for Ce<sup>3+</sup>/ $Tb^{3+}$ -doped phosphor<sup>5,6</sup>.

Scattering, direct and indirect excitation emission spectra of a Eu<sup>2+</sup>-doped thin film and PbS quantum dots can be seen in Figure 3 and 4, respectively. Note the considerable emission measured under indirect excitation of the thin film in the integrating sphere.







Figure 2: Scattering and emission spectra of Ce<sup>3+</sup>/Tb<sup>3+</sup>-doped phosphor. The measurement conditions were:  $\Delta\lambda_{exc}=10$  nm,  $\Delta\lambda_{em}=0.5$  nm, step=0.25 nm, dwell=0.1s.

## ABSOLUTE QUANTUM YIELD OF UV-TO-NIR EMITTING SAMPLES



Figure 3: Eu<sup>2+</sup>-doped thin film under (a) direct and (b) indirect excitation. The samples were excited at 480 nm with excitation and emission bandwidths  $\Delta \lambda_{exc}$ =3 nm and  $\Delta \lambda_{em}$ =0.2 nm, respectively, with a step of 0.25 nm and integration time of 0.1s.



Figure 4: Scattering and emission spectra of PbS quantum dots. Conditions were:  $\Delta \lambda_{exc} = 5 \text{ nm}, \Delta \lambda_{em} = 2 \text{ nm}, \text{ step} = 1 \text{ nm}, \text{ dwell} = 0.5 \text{ s}.$ 

#### REFERENCES

- Crosby, Glenn A., and James N. Demas. "Measurement of photoluminescence quantum yields. Review." The Journal of Physical Chemistry 75.8 (1971): 991-1024.
- Würth, C., Grabolle, M., Pauli, J., Spieles, M. & Resch-Genger, U. Comparison of Methods and Achievable Uncertainties for the Relative and Absolute Measurement of Photoluminescence Quantum Yields. Anal. Chem. 83, 3431–3439 (2011).
- Würth, C., Geißler, D., Behnke, T., Kaiser, M. & Resch-Genger, U. Critical review of the determination of photoluminescence quantum yields of luminescent reporters. Anal Bioanal Chem 407, 59–78 (2014).
- de Mello, J. C., Wittmann, H. F. & Friend, R. H. An improved experimental determination of external photoluminescence quantum efficiency. Adv. Mater. 9, 230–232 (1997).
- Jüstel, T., Krupa, J.-C. & Wiechert, D. U. VUV spectroscopy of luminescent materials for plasma display panels and Xe discharge lamps. Journal of Luminescence 93, 179–189 (2001).
- Bühler, G. & Feldmann, C. Microwave-Assisted Synthesis of Luminescent LaPO<sub>4</sub>:Ce,Tb Nanocrystals in Ionic Liquids. Angewandte Chemie International Edition 45, 4864–4867 (2006).

#### For more information, contact:

#### Edinburgh Instruments

2 Bain Square, Kirkton Campus, Livingston, EH54 7DQ United Kingdom.

T: +44 (0)1506 425 300 F: +44 (0)1506 425 320

E: sales@edinst.com

W: www.edinst.com



#### +44 (0)1506 425 300 | sales@edinst.com | www.edinst.com

Words and logos marked with ® or ™ are registered trademarks or trademarks owned by EDINBURGH INSTRUMENTS Limited. Other brands and names mentioned herein may be the trademarks of their respective owners. Neither the whole nor any part of the information contained in, or the product described in, this document may be adapted or reproduced in any material form except with the prior written permission of the copyright holder.