

## PHOTO-PHYSICAL PROPERTIES OF ANCIENT AND MODERN ARTWORK PIGMENTS

AN\_P23 v.2; Dec. 15, Gianluca Accorsi, Giovanni Verri, Georgios Arnaoutakis

EDINBURGH  
INSTRUMENTS

## INTRODUCTION

The interest in non-invasive investigations for the study and identification of painting materials has greatly increased during the past two decades.<sup>1</sup> The nature of museum objects is such that sampling is always kept to a minimum, thus, a spectroscopic approach is ideal for this purpose. In this respect, an ancient and a modern pigment, cuprorivaite,



FLS980 Fluorescence Spectrometer

$\text{CaCuSi}_4\text{O}_{10}$  (Egyptian blue)<sup>2</sup> and barium manganate (VI) sulphate compound (manganese blue)<sup>3</sup> respectively, have been photo-physically investigated. The investigation includes qualitative information, i.e. emission and excitation spectra, as well as quantitative that is the photoluminescence quantum yields (PLQY) that have been obtained in the near-infrared (NIR) region.

## METHODS &amp; MATERIALS

Excitation and emission spectra were measured using an FLS980 Fluorescence Spectrometer equipped with a 450 W Xe lamp with double excitation and emission monochromators. Near-infrared detectors (Hamamatsu) were used for the detection of both samples. PLQY have been calculated by corrected emission spectra obtained by using a barium sulphate coated integrating sphere, following the procedure described by De Mello et al.<sup>4</sup> Experimental uncertainties were estimated to be  $\pm 20\%$  for emission quantum yields,  $\pm 2$  nm

## RESULTS - DISCUSSION

and  $\pm 5$  nm for absorption and emission peaks, respectively.

Figure 1 reports the excitation and luminescence spectra of the two samples. Egyptian blue shows two different electronic transitions ( ${}^2B_{1g} \rightarrow {}^2E_g$  and  ${}^2B_{1g} \rightarrow {}^2A_{1g}$ ) that can be assigned to  $\text{Cu}^{2+}$  ions, which are expected to be the only luminescent components of cuprorivaite.<sup>1</sup> On the other hand, those of ManganeseBlue can be attributed to ligand-field (LF) transitions ( ${}^2E \rightarrow {}^2T_2$ , c.a. 800 nm - 900 nm) and to ligand-to-metal charge-transfer bands (between 600 nm and 800 nm) of the  $\text{MnO}_4^{2-}$  unit.<sup>5</sup>

The luminescence profiles have peaks at 920 nm and 1300 nm for Egyptian blue and manganese blue, respectively.

Consequently, the corresponding PLQYs have been obtained.

Manganese blue showed a quantum yield of  $\Phi = 0.5\%$ <sup>6</sup> attributed to luminescence quenching factors that occur at lower energy. By contrast, Egyptian blue revealed to be a very strong NIR emitter ( $\Phi = 10.5\%$ )<sup>1</sup> that, to the best of our knowledge, has the highest quantum efficiency for a molecule-level chromophore in the 800 nm – 1100 nm range.

## CONCLUSION

The photo-physical properties of an ancient and a modern art pigment were investigated by means of fluorescence spectroscopy. In addition, the PLQY of the pigments has been measured, as a quantitative means of their photoluminescence properties and demonstrate the effectiveness of this method as a minimally invasive means of measuring samples of cultural significance.




Figure 1: Excitation and emission spectra of ancient, Egyptian blue, and modern, manganese blue, pigments. The energy levels corresponding to each transition are also displayed.



FLS Spectrometer Integrating Spheres

## REFERENCES



1. Accorsi, G.; Verri, G.; Bolognesi, M.; Armaroli, N.; Clementi, C.; Miliani, C.; Romani, A. The exceptional near-infrared luminescence properties of cuprorivaite (Egyptian blue). *Chem. Commun.* 2009, (23), 3392-3394.
2. (a) *Artist's Pigments – A Handbook of Their History and Characteristics*. Elisabeth West Fitzhugh editor: 1997; Vol. 3; (b) Berke, H. The invention of blue and purple pigments in ancient times. *Chem. Soc. Rev.* 2007, 36 (1), 15-30.
3. Eastaugh, N.; Walsh, V.; Chaplin, T.; Siddall, R. *Barium Manganate(IV) Sulfate*. Butterworth-Heinemann: Oxford (UK), 2008.
4. De Mello, J. C.; Wittmann, H. F.; Friend, R. H. An improved experimental determination of external photoluminescence quantum efficiency. *Adv. Mater.* 1997, 9 (3), 230-236.
5. Chen, Z. L.; Schmale, H. W.; Fox, T.; Blacque, O.; Berke, H. NH-functionalized tungsten complexes of 2-(dimethylphosphino)imidazole. *J. Organomet. Chem.* 2007, 692 (22), 4875-4885.
6. Accorsi, G.; Verri, G.; Acocella, A.; Zerbetto, F.; Lerario, G.; Gigli, G.; Saunders, D.; Billinge, R. Imaging, photophysical properties and DFT calculations of manganese blue (barium manganate(VI) sulphate) - a modern pigment. *Chem. Commun.* 2014, 50 (97), 15297-15300.

## For more information, contact:

## Edinburgh Instruments

2 Bain Square,  
Kirkton Campus,  
Livingston, EH54 7DQ  
United Kingdom.

T: +44 (0)1506 425 300

F: +44 (0)1506 425 320

E: [sales@edinst.com](mailto:sales@edinst.com)

W: [www.edinst.com](http://www.edinst.com)

