LP980

Transient Absorption Spectrometer





The LP980 spectrometer sets the standard for technical performance required in a transient absorption research instrument, offering unsurpassed measurement capabilities across a broad range of applications.

A first of its kind in the market, the LP980 allows for the measurement of transient absorption and both laser-induced fluorescence and phosphorescence up to 2.55 µm, from nanoseconds to seconds. Spectral and kinetic acquisition are possible in the same instrument thanks to a dual monochromator/spectrograph.

Measurements:

- > Transient Absorption
- > Laser-Induced Fluorescence (LIF) and Phosphorescecene
- > Time-Gated Spectra (ns, ms, ms time ranges)
- > Triplet-Triplet Annihilation
- > Reaction Rate Studies

Key Features



Dual sample chamber

Measurement of transient absorption and laser-induced fluorescence



Two detectors

Spectral (ICCD) and/or time-resolved (PMT) data



ns to seconds

Accurate kinetics over 9 orders of magnitude



High sensitivity

Detection limit of Δ OD 0.002 (kinetic) and Δ OD 0.0005 (spectral)



L900 software

One comprehensive package for control of all components and measurements



LP980 - Base Configuration

The LP980 is a transient absorption spectrometer using the pump-probe technique for measuring transient kinetics (Kinetic Mode) and/or time-gated transient spectra (Spectral Mode), generated by laser excitation.

Transverse sample excitation geometry comes as standard. Thin-film, diffuse reflection, fluorescence and phosphorescence lifetime measurement, accessories are available as options.

MONOCHROMATOR / SPECTROGRAPH	Туре	Czerny-Turner with triple grating turret
	Focal length	325 mm
	Mirror	Automatic, computer-controlled for detector selection
	Slits	5 mm to 10 mm (continuously adjustable), motorised
LASER EXCITATION SOURCE*	Single wavelength	Flashlamp pumped Q-switched Nd:Yag laser operating at 1064 nm, 532 nm, 355 nm, or 266 nm*
	Tuneable	OPO, tuneable in range 410 nm – 710 nm (signal). Idler and UV doubler options possible

^{*} We can supply a fully integrated laser, please contact us for more information.

LP980-K (Kinetic Mode)

For lifetime transient decay measurements at a single wavelength.

GRATING Plane ruled grating 1800 grooves/mm, 500 nm blaze as standard

 DISPERSION
 1.66 nm/mm

 SPECTRAL RANGE
 200 nm - 870 nm

SPECTRAL RESOLUTION 0.1 nm

SENSITIVITY Δ OD 0.002 (single shot - fast detector option, PMT), Δ OD 0.0005 (single shot - slow detector option, ICCD)

DETECTOR TYPE Photomultiplier with 5 stage dynode chain for high current linearity

DETECTOR IMPEDANCE 50 Ω (amplified – fast detector, <3 ns rise time), 1 k Ω (slow detector, <100 μ s rise time)

LP980-KS (Kinetic & Spectral Mode)

For lifetime transient decay measurements AND spectral measurements of the decay process.

GRATING Kinetic Mode grating plus an additional plane ruled grating: 150 grooves/mm, 500 nm blaze supplied

DISPERSION 19.9 nm/mm

SPECTRAL COVERAGE 520 nm (active horizontal ICCD dimension: 25 mm)

SPECTRAL RESOLUTION 0.52 nm (spectral coverage / 960 pixels)

SENSITIVITY Δ OD 0.0005 (single shot)

DETECTORSKinetic Mode PMT plus an additional image intensified CCD camera (ICCD) supplied

 MIN. OPTICAL GATE WIDTH
 7 ns (FWHM)

 ACTIVE PIXELS
 960 x 256

 ACTIVE AREA
 25 mm x 6.7 mm

COOLING -20°C as standard (-30°C with additional water circulation)

LP980 Upgrade Options

GRATING OPTIONS A variety of gratings are available with 150-2400 grooves/mm, optimised from UV through to NIR

SAMPLE HOLDER OPTIONS Cross-beam geometry, diffuse reflectance, thin-film, LIF

DETECTOR OPTIONS PMT-980 (200nm - 980 nm), InGaAs Detectors (900 nm - 2550 nm), NIR-PMT (up to 1650 nm)



