APPLICATION NOTE

Batch Processing and Global Analysis of Fluorescence Lifetimes

AN_P83; Georgios Arvanitakis

Introduction

Analysing large fluorescence lifetime datasets requires automated processing and the identification of common patterns across the dataset. Batch Processing and Global Analysis are two automated processing methods for fluorescence lifetime datasets that are available in Edinburgh Instruments FAST lifetime analysis software. In this Application Note, fluorescence lifetimes and relative weights of a threefluorophore mixture of erythrosin b, fluorescein, and 9aminoacridine (9AA) were analysed using both batch processing and global analysis in FAST, and the processing methods compared.

Methods and Experimental Setup

The fluorescence lifetimes were acquired using an Edinburgh Instruments FLS1000 Photoluminescence Spectrometer equipped with a PR2 Microwell Plate Reader and a High-Speed Hybrid Photodetector (HS-HPD). An EPL-445 picosecond pulsed diode laser was used as the excitation source and operated at 10 MHz, and the lifetimes were acquired using TCSPC. Each fluorophore was diluted in methanol to give an absorbance of around 0.03 for erythrosin b, 0.06 for fluorescein, and 0.09 for 9AA at the EPL's operational wavelength, 437 nm. Three wells were filled with ~150 μ l of each pure fluorophore mixture. Representative figures of the plate reader and the plate reader interface in Fluoracle that shows the ten selected wells containing the mixture can be seen in Figures 1a and 1b, respectively.




Figure 1 a) The FLS1000 plate reader, and b) PR2 Plate Reader Interface in Fluoracle.

Lifetime Analysis of the Individual Decays

The fluorescence lifetimes of the three pure fluorophores can be seen in Figure 2. Erythrosin b (Figure 2a) exhibits the shortest decay (0.46 ns), and therefore, the fit was a deconvolution of the



instrument response function (IRF) and the exponential decay of the sample. The IRF was not considered for Fluorescein and 9AA as they were characterised by longer lifetimes. Figures 2b and 2c show the Fluorescein (4.41 ns) and 9AA's (12.41 ns) decays, fitted by a single exponential decay.




Figure 2 Fluorescence lifetimes of a) erythrosin b, b) fluorescein, and c) 9AA.

Batch Processing Analysis of the Mixture

The ten fluorescence decays of the mixture were first analysed using batch processing. In batch processing, multiple decays can be automatically analysed in a single batch. The decays are analysed independently and fitted to a chosen model to obtain independent information for each decay. The ten decays of the mixture were imported into FAST, and batch processed using the Exponential Components Analysis routine. A threecomponent exponential model was chosen since the mixture comprises three fluorophores (Figure 3). In the case of an unknown sample, the user could start fitting with a single exponential component and then add additional components if required after reviewing the residuals and chi-squared optimisation. FAST then performs fitting for each fluorescence decay providing numerical information about τ_1 , τ_2 , and τ_3 , their amplitudes, B_i , (with amplitude standard deviation, ΔB_i) and weights, f_i (%) (with weights standard deviations, Δf_i %).

APPLICATION NOTE

Batch Processing and Global Analysis of Fluorescence Lifetimes



AN_P83; Georgios Arvanitakis

Tail fitting	O Reconvoluti	on
Automated fitting	🖉 Manual para	am. contro
From 0	To 500	[Chnls
For fit from peak	set "From" equa	10.
Number of Compo	nents 3 🔒	
Amplitude Fi	x Lifetime	Fix
B ₁ 3333.3	τ ₁ 0.4254	
B ₂ 3333.3	τ ₂ 1.3452	
B ₃ 3333.3	τ ₃ 4.2540	
B4] T4	
	Value	Fix
Shift	(ns) 0.000	
Decay backgr	ound 0.6300	
IRF backgr	ound	
Supp	ort Plane Analysis	s 🗆

Figure 3 Parameters for batch processing in FAST.

The distributions of the relative weights, f_i (%), of each exponential component are plotted in Figure 4 alongside the average lifetime of each component. The average lifetimes of the τ_1 , τ_2 and τ_3 components closely match those of pure erythrosin B, fluorescein and 9AA, respectively.




Figure 4 Distribution of the lifetimes and relative weights of the three components across the ten decays using batch processing.

Global Analysis of the Mixture

In global analysis, the component lifetimes are no longer independent across the multiple decays in the dataset; they are

a global parameter. Global analysis assumes that a global component lifetime can describe all the fluorescence decays in a dataset.

Tail fitting	O Reconvolution
Standard Global F	it 📿 Advanced Global
From 0	To 500 [Chnl
For fit from pea	k set "From" equal 0.
Number of Compo	inents 3 🕂
Amplitude Fix L	ink Lifetime Fix L
B1 3.33e+0(τ ₁ 0.425
B2 3.33e+0(τ ₂ 1.35
B3 3.33e+0(τ ₃ 4.25
B4 🗌 🗆 [⊐ T4 🔽 🗆 [
	Value Fix
Shift	(ns) 0.000
Decay Backg	round 0.63
IRF Backg	round

Figure 5 Parameters for global analysis in FAST.

Figure 5 shows how to set the parameters in FAST to perform global analysis. τ_1 , τ_2 and τ_3 are set to be linked. Linking a lifetime component means that the same component lifetime is a global parameter and is used to fit all decays in the dataset. The component's lifetime is then optimised across all decays in the dataset during the fitting routine to find the lifetime that best describes all decays collectively.




Figure 6 Lifetimes and relative weights of the three components using global analysis across the ten decays.

APPLICATION NOTE

Batch Processing and Global Analysis of Fluorescence Lifetimes

AN_P83; Georgios Arvanitakis



The lifetimes and relative weight of the global analysis results are shown graphically in Figure 6. Since the lifetimes are now a global parameter, there is no longer a distribution for each lifetime component. Instead, there is a single lifetime value for each component that best describes all ten decays. It is worth mentioning that if τ_1 , τ_2 , and τ_3 were all set to unlinked, then the results would precisely match the batch processing analysis of the previous section. Thus, the user can perform both techniques promptly by linking and unlinking the time variables.

Conclusion

This Application Note outlined the successful application of Edinburgh Instruments' FAST software for batch processing and global analysis in studying the lifetimes of a three-fluorophore mixture. The batch analysis approach treats the lifetimes of each component in the mixture as independent parameters, while the global analysis approach considers them as linked (global parameters). Both methods offer valuable insights into the lifetimes and relative weight distribution of the individual components within a multi-fluorophore mixture. The choice between these methods should be based on the specific samples being studied and their intended applications.



For more information, please contact:

+44 (0) 1506 425 300 sales@edinst.com www.edinst.com