


This mini remote controlled, single point magnetic stirrer has been designed such that it easily fits underneath the standard Edinburgh Instruments single cuvette holder. The stirrer comprises three stirrer bars and a stand alone stirrer controller. A long, highly flexible 0.3mm thick foil cable joins the magnetic stirrer to the separate remote control unit. The remote control nature of this device and the large distance between the magnetic stirrer and the stirring vessel means that it is ideally suited for integration into our analytical spectrometers.



The unit features the following:

- Wear free inductive VARIOMAG drive
- Remote controlled by separate control unit
- No larger than 1 Euro coin
- Ideally suited for integration into the standard Edinburgh Instruments single cuvette holder
- Highly flexible foil cable, 0.3mm thick
- Immersible, can be used directly in the water bath
- Safe low voltage



Edinburgh Instruments has a policy of continuing product development and reserve the right to amend specifications without prior notice. Dec. 2005